



Z-Band[®] UltraSoft Wristbands



SEE MORE. DO MORE.



Healthcare providers can take a softer approach to patient identification with Zebra's patent-pending Z-Band[®] UltraSoft direct thermal adhesive wristband.

Featuring a soft, flexible material, Z-Band UltraSoft is the softest direct thermal wristband in the market. Z-Band UltraSoft also offers excellent resistance to most hand sanitizers, supporting patient safety by ensuring readable text and scannable barcodes for up to 14 days.

Zebra's Z-Band UltraSoft direct thermal wristbands allow you to print patient information and barcodes directly on the wristbands.

Ideal for These Applications

Patient Identification / Hospitals, Outpatient clinics

A barcode that fails to scan or requires repeated attempts to scan jeopardizes patient safety and staff productivity. Zebra's supplies research team tested the image durability of Z-Band direct thermal wristbands to ensure optimum image durability and scanning performance.





For more details on the durability laboratory testing performed, please download our white paper on www.zebra.com titled "*Laboratory Testing of Zebra's Z-Band[®] Direct Thermal Antimicrobial Wristbands.*"

THE SOFTEST DIRECT THERMAL WRISTBAND IN THE MARKET

Z-Band UltraSoft:

- Ultra-Soft wristband.
- Resistance to hand sanitizers such as Purell®.
- Adhesive closure.
- Scannable up to 14 days.

Z-Band UltraSoft is available in cartridges for the HC100™ wristband printer in the following sizes:

PART NUMBER	SIZE	
ADULT		
10015355K	1" x 11"	
10015356K	0.75" x 11"	
PEDIATRIC		
10015357K	1" x 7"	
INFANT		
10015358K	1" x 6"	

In addition, Zebra's Z-Band direct thermal wristbands:

- Feature a silver antimicrobial coating that protects the wristband.
- Are scannable for up to 14 days.
- Are quick and easy for staff to print and securely fasten in less than a minute using an adhesive tab or color clip.
- Are made of polypropylene and vinyl. The direct thermal imaging chemistry is composed primarily of leuco dyes and developers. The inks/varnishes used are UV varnishes and water based inks. **There is no metal content in or on these particular patient identification bands.** Therefore, Z-Band UltraSoft wristbands are considered to be MR-Safe.



Corporate Headquarters
+1 800 423 0442
inquiry4@zebra.com

Asia-Pacific Headquarters
+65 6858 0722
apacchannelmarketing@zebra.com

EMEA Headquarters
+44 (0)1628 556000
mseurope@zebra.com

Latin America Headquarters
+1 847 955 2283
inquiry4@zebra.com

Other Locations / USA: California, Georgia, Illinois, Rhode Island, Texas, Wisconsin **Europe:** France, Germany, Italy, the Netherlands, Poland, Spain, Sweden, Turkey, United Kingdom **Asia Pacific:** Australia, China, Hong Kong, India, Indonesia, Japan, Malaysia, Philippines, Singapore, South Korea, Taiwan, Thailand, Vietnam **Latin America:** Argentina, Brazil, Colombia, Florida (LA Headquarters in USA), Mexico **Africa/Middle East:** Dubai, South Africa